

RNA Extraction

Reliable RNA purification

Total RNA Extraction Kit (Blood/Bacteria/Cultured Cells)

The Total RNA Extraction Kit (Blood/Bacterial/Cultured Cells) is specially designed for purification of total RNA from bacterial, cultured cells and fresh human whole blood. The method uses detergents and a chaotropic salt to lyse the cells and inactivate RNase. The lysate is applied to provided Lysate Filter columns (Mini System) to remove cell debris and ensure complete cell lysis. RNA in chaotropic salt solutions binds to the glass fiber matrix of the columns. Following washing off of contaminants, the purified RNA is eluted by RNase-free water. ssRNA and dsRNA of > 200 bps to 1000's of bps in length are efficiently purified.

Applications


RT-PCR, cDNA synthesis, In Vitro Translation, Northern blotting, Micorarray prep

Quality Control

The quality of Total RNA Extraction Kit (Blood/Bacterial/Cultured Cells) is tested on a lot-to-lot basis. The kits are tested by isolation of total RNA from 300µl of fresh human whole blood. More than 1µg of total RNA was quantified with a spectrophotometer and checked by formaldehyde agarose gel analysis. Finally, RT-PCR was used to ensure the quality of total RNA.

Protocol

Sample (Protocols Included*):	Animal Tissues, Whole Blood/ Buffy Coat , Cultured Animal Cells, Gram -ve/+ve Bacteria
Genomic RNA Extraction Kits <i>Mini</i>	
Sample size:	106 Animal Cells, 300µl Whole Blood, 108 Bacterial Cells
Preparation time:	20 min
Yield:	< 30 µg
Genomic RNA Extraction Kits <i>Maxi</i>	
Sample size:	100-200mg Animal Tissue, 107-108 Cultured Cell, 5ml Blood Sample
Preparation time:	60 min
Yield:	< 500 µg
*For certain protocols items such as DNase I, Lysozyme may need to be purchased separately.	



Storage and Stability

This kit should be stored at room temperature.

Reference

Vogelstein, B., and Gillespie, D. (1979) Proc. Natl. Acad. Sci. USA 76, 615.

Additional Requirements(dependent on sample type):

Blood Protocol:

PBS, 70% Ethanol and 96-100% Ethanol, Sterile, RNase-free pipet tips and microcentrifuge tube DNase I: DNase I (RNase free), RNase-free -Mercaptoethanol (b-ME) must be added to RB Buffer before use.

Bacterial Protocol:

Lysozyme Buffer: 20 mg/ml lysozyme; 20 mM Tris-HCl; 2 mM EDTA; 1% Triton X-100; pH 8, Prepare the lysozyme buffer immediately prior to use.

Cat.No. YRB50

50 mini preps / kit
RBC Lysis Buffer: 120 ml
RB Buffer: 30ml
RT Buffer: 15 ml
R-W1 Buffer: 25 ml
R-Wash Buffer (concentrated): 25 ml
RNase-Free Water: 10ml
Filter Column Set: 50 sets
(Comes with 2ml Collection Tube)
RB Column Set: 50 sets
(Comes with 2ml Collection Tube)


Cat.No. YRB100

100 mini preps / kit
RBC Lysis Buffer: 120ml x2
RB Buffer: 60ml
RT Buffer: 30ml
R-W1 Buffer: 50ml
R-Wash Buffer (concentrated): 25ml
RNase-Free Water: 10ml
Filter Column Set: 100 sets
(Comes with 2ml Collection Tube)
RB Column Set: 100 sets
(Comes with 2ml Collection Tube)

All components are RNase-Free treated
RB Column Set is blister packed

All components are RNase-Free treated

Performance



RBC Bioscience Labs

1: mouse liver (6µg)

2: kidney (6µg)

3: cos7 cells (6µg)

4: hela cells (6µg)



Quality Nucleic Acid Purification Systems

RBC Real Genomics DNA/RNA Purification



RBC Bioscience Corp.
www.rbcbioscience.com
info@rbcbioscience.com