



Important Notes For MagCore[®] triXact RNA Kit



Please subject fresh samples (whole blood, cultured cells and tissues) to the protocol **within 24 hours from obtaining them**. Or immediately lyse the samples in RB buffer for stabilization, then store them at -80°C for up to 1 month.



During sample preparation, **keep the sample on ice** for as long as possible to avoid RNA degradation.



When doing extraction with MagCore[®] triXact RNA Kit, always wear a suitable lab coat, disposable gloves, and protective mask. Do not talk during the experiment to avoid contamination.



Ensure that the experimental environment is suitable for operating RNA experiments. Performing the extraction in the hood is recommended.