



## Lysozyme Set

Cat.No.YLYS250  
1000mg  
Activity:>23.5 KU/mg

V1.0

### Kit contents

Content	
Lysozyme	250mg/vial *4
Lysozyme Reaction Buffer	60 ml*1

### Shipping and Storage

The Lysozyme Set is shipped at room temperature and should be stored immediately upon receipt at 2–8 °C up to 12 months.



### Unit Definition

One unit produces a decrease the absorbance at 450 nm by 0.001 per minute at pH 7.0 and 25°C.  
Use a 0.3 mg/ml suspension of *Micrococcus lysodeikticus* cells as a substrate.

### Preparation before using

Freshly prepared 20 mg/ml Lysozyme solution with Lysozyme Reaction Buffer before use.

### Application

Nucleic Acid Extraction, Detection of BSE forming proteins which are uniquely resistant to proteolytic cleavage  
Analysis of membrane structure.



[www.rbcbioscience.com](http://www.rbcbioscience.com)



RBCBioscience